



Isolation of quartz for cosmogenic in situ ¹⁴C analysis

Keir A. Nichols¹, Brent M. Goehring¹

¹Department of Earth and Environmental Sciences, Tulane University, New Orleans, LA, 70118, USA *Correspondence to*: Keir A. Nichols (knichol3@tulane.edu)

- 5 Abstract Froth flotation is a commonly used procedure for separating feldspars and micas from quartz for the preparation of quartz mineral separates to carry out cosmogenic nuclide analysis. Whilst extracting carbon from quartz we observed in situ carbon-14 (¹⁴C) concentrations which were anomalously high and in excess of theoretical geological maximum concentrations. Further etching of sample material reduced carbon yields and ¹⁴C concentrations, yet the latter remained unrealistically high. When quartz from the original whole rock sample was isolated in our laboratory, we observed even lower carbon yields and
- 10 geologically plausible in situ ¹⁴C concentrations. After ruling out unlikely geological scenarios and systematic measurement issues, we decided to investigate the quartz isolation procedure as a potential source of ¹⁴C contamination. We hypothesised that laurylamine (dodecylamine), an organic compound used as part of the froth flotation procedure, elevates ¹⁴C concentrations if residual laurylamine is present. We demonstrate that laurylamine has a ¹⁴C modern carbon source and thus has the potential to influence in situ ¹⁴C measurements if present in minute but measurable quantities. Furthermore, we show
- 15 that insufficient sample etching results in laurylamine-derived carbon persisting through step heating of quartz and is subsequently collected with the in situ component released at 1100 °C. We therefore demonstrate that laurylamine contaminates in situ ¹⁴C measurements. We provide guidelines for the preparation of quartz based on methods developed in our laboratory and demonstrate that all laurylamine derived carbon and ¹⁴C is removed when applied. We recommend that the procedures presented be used at a minimum when using froth flotation to isolate quartz for in situ ¹⁴C measurements.

20 1.0 Introduction

In the course of extracting carbon from quartz we have, on multiple occasions, observed concentrations of in situ ¹⁴C that were anomalously high and in excess of geologically plausible maximum concentrations. We hypothesise that the elevated in situ ¹⁴C concentrations are sourced from part of the widely used mineral separation procedure known as froth flotation, a process that relies on the long-chain organic compound laurylamine. Our observations, combined with a desire to continue use

25 of froth flotation for the benefits it provides during quartz separation, form the motivation for this paper. In this study we explore both the potential influence that quartz isolation procedures have on resulting measured ¹⁴C concentrations as well as procedures to minimise potential contamination during use of froth flotation.

Froth flotation is a method by which feldspars and, to a lesser extent, micas are separated from quartz (Herber, 1969). The method precedes acid etching as part of the quartz isolation process for cosmogenic nuclide analysis and is used by





- 30 numerous cosmogenic nuclide laboratories. It is useful for samples containing large proportions of feldspars. The method vastly reduces the resources required to etch samples, and also helps avoid the relatively laborious and time-consuming method of density separation. A motivating factor for this study was the realisation that froth flotation introduces carbon to sample material through the use of the organic compound laurylamine (also known as dodecylamine, C₁₂H₂₇N). This was previously of no consequence as the method was primarily used to isolate quartz for the measurement of ¹⁰Be and ²⁶Al. There is no
- 35 standard procedure for froth flotation or post-froth flotation sample etching. Different laboratories use various quantities of laurylamine and varying etching procedures, which complicates the matter further. The carbon content, and especially the ¹⁴C content, of laurylamine has yet to be measured, thus the potential for contamination of in situ ¹⁴C measurements is unquantified. In the first part of this work, we summarise the froth flotation procedure as well as the overall quartz isolation process
- used for in situ ¹⁴C analysis. We describe the range of methodologies used today, and detail those used at Tulane. We then describe the initial measurements that led us to hypothesise that froth flotation could be causing contamination of in situ ¹⁴C results. Finally, we describe the methodology and results of a systematic study that demonstrates laurylamine contains modern carbon, that froth flotation does contaminate samples with regards to both ¹⁴C and carbon in general, and that contaminant ¹⁴C can be removed with sufficient sample etching. We demonstrate that the post-froth flotation etching methodology used in our laboratory ensures that quartz is isolated effectively and without influencing the resulting in situ ¹⁴C measurements. We conclude that froth flotation should be applied with care if in situ ¹⁴C is to be measured, and that the post-froth flotation etching
- methodology described below should be applied at a minimum to ensure that samples are free of contaminant ¹⁴C from laurylamine.

1.1 Froth flotation and the isolation of quartz from whole rock material

50 1.1.1 Pre-Froth Flotation

55

Prior to froth flotation, whole rock material is typically crushed, milled and sieved to isolate the 250 - 500 µm size fraction. This is then rinsed with tap or deionized water to remove any fine grain-sized material. At this point samples are ready for froth flotation, although we commonly first dry samples so that a magnetic mineral separation can be performed to remove any mafic material present prior to frothing, which we find improves overall frothing efficiency. The sample is ready for froth flotation following the removal of fine grain-sized material and the optional magnetic separation.

method Our for froth flotation is largely based on that used at PRIME Lab (http://www.physics.purdue.edu/primelab/MSL/froth floatation.html). The first stage of froth floatation is the conditioning of sample material with dilute (< 5% v/v) hydrofluoric acid (HF). Conditioning the sample makes the feldspar (and mica) grains hydrophobic and the quartz grains hydrophilic, which is key to the separation process. We condition each sample in a 1 L

60 Nalgene bottle with enough 5 % HF/HNO₃ to saturate and cover the sample, without agitation beyond gently swirling the bottle a few times. The sample is left to sit for no more than five minutes before decanting the acid solution and beginning





froth flotation. Some laboratories condition the sample with dilute HF (1 to 5 %) for up to 60 minutes on a shaker table; we have found that five minutes with 5 % HF/HNO₃ is sufficient and improves separation efficiency.

65 1.1.2 Frothing Solution

Laurylamine is combined with glacial acetic acid by dissolution to form the frothing solution, typically at a 1:1 ratio. We combine approximately 300 ml of glacial acetic acid and 300 ml of laurylamine to form a 600 ml stock frothing solution. The frothing solution is then combined with water and carbonated or mixed with bubbly tap water. Some laboratories add the concentrated frothing solution directly to sample material, followed by the addition of carbonated or bubbly tap water. Other 70 laboratories make the stock solution of glacial acetic acid and laurylamine and combine it with water before adding it to the sample material. For each sample, we combine approximately 6 ml of frothing solution with 20 L of water. The net concentration of both acetic acid and laurylamine in the frothing solution is 0.03% v/v. In terms of their purpose in the froth flotation procedure, laurylamine acts as a collector agent, or surfactant, and is thus required to separate the hydrophobic and hydrophilic mineral grains. Glacial acetic acid is used because laurylamine dissolves into it more readily compared to water, and it keeps the pH of the solution low.

75

1.1.3 Froth Flotation

A few drops of eucalyptus oil are added to the sample in a metal bowl before the frothing solution is applied to the sample. The eucalyptus oil holds the bubbles together to which the feldspar and mica grains attach. We use a hose connected to a soda-fountain carbonator to dispense the frothing solution. The now carbonated and dilute frothing solution is used to 80 move the sample material from the 1 L bottle to the metal bowl. The frothing solution is then applied to the sample material in the bowl. The feldspar grains, owing to their hydrophobic nature, float to the top of the mixture whilst the quartz grains remain at the bottom. We apply 3 to 4 L of dilute frothing solution to the sample before waiting a few seconds and decanting the feldspar grains into a second metal bowl. The feldspar grains are usually discarded, though they may be saved for ³⁶Cl analysis. The froth flotation procedure is repeated until most of the feldspar fraction has been removed or no additional 85 separation of quartz and feldspar is accomplished. For a granitic sample of ca. 400 to 500 g, we find that five to six rounds of

froth flotation are needed before either the froth flotation process is complete and the vast majority of feldspar has been removed, or froth flotation becomes less effective and the sample requires further conditioning. After reconditioning the sample in 5 % HF/HNO₃ for five minutes, additional rounds of froth flotation can be performed.

90





1.1.4 Post-Froth Flotation Acid Etching

Froth flotation is followed by etching the sample in HF or HF/HNO₃ to remove extraneous minerals, to partially dissolve or etch the quartz grains to remove meteoric cosmogenic nuclides, and to lower major ion concentrations (e.g., Fe, Ti, Al). Generally, the etching process follows a heavily modified version of the method of Kohl and Nishiizumi (1992). A 95 typical procedure used by many laboratories involves first etching samples in 5 or 1 % HF or HF/HNO₃ on a shaker table or sample roller for multiple periods, followed by etching in 1 % HF or HF/HNO3 in an ultrasonic bath. Some laboratories etch samples in an ultrasonic bath without the use of a sample roller or shaker table beforehand. Between etches, samples are rinsed with deionised water (i.e., 18.2 MQ H₂O). The number of etches will vary due to factors such as sample lithology, amount of sample material, effectiveness of the froth flotation procedure, in addition to the varying minimum standard procedures for a given laboratory. For in situ ¹⁴C analysis, samples are often etched until they pass a visual test under a binocular microscope 100 and the sample appears to be solely composed of quartz.

1.2 Initial Anomalous C-14 Measurements

Whilst measuring the in situ ¹⁴C concentration of glacial erratic samples as part of multiple projects we observed in situ ¹⁴C measurements that were in excess of geologically plausible maximum concentrations (Fig. 1). In each case the 105 maximum concentration for a sample is set by the in situ ¹⁴C saturation concentration for the given sample location, shown in Fig. 1. The only way elevated in situ¹⁴C concentrations could be explained is with an unlikely geomorphic scenario in which the samples were exposed at much higher elevations for a significant period of time before being rapidly transported to their sampling location. This phenomenon was described by Balco et al. (2016) and potentially observed by Balco et al. (2019). Subsequent elevated in situ ¹⁴C concentrations measured from bedrock samples led us to rule out this scenario as the sole 110

source of the observed elevated in situ ¹⁴C concentrations, and we began to explore other explanations.

To investigate the cause or causes for the anomalously high in situ ¹⁴C measurements we performed additional measurements of in situ¹⁴C concentrations of samples displaying elevated concentrations following additional etches in 1 % HF/HNO₃ for two 24-hour periods. This yielded lower, but still anomalously high, ¹⁴C concentrations (Fig. 1), as well as lower

- total carbon yields (Hillenbrand, unpub.). We note that the quartz from which the anomalously high in situ¹⁴C concentrations 115 and elevated carbon yields were measured was isolated at other laboratories that use slight variations in their quartz isolation procedures to ours. To investigate further, we measured the in situ ¹⁴C concentration from the same samples but isolated the quartz from whole rock material using our standard procedure (Sect. 1.1). Carbon yields were further reduced, and the resulting ¹⁴C concentrations were both lower and geologically plausible (Fig. 1).
- 120

The additional measurements left two potential explanations for the elevated concentrations; unidentified systematic measurement issues or contamination of sample material. Repeat measurements of the quartz interlaboratory comparison material CRONUS-A (Jull et al., 2015; Goehring et al., 2019) and other samples allowed us to rule out systematic measurement





125

procedure was a potential source of ¹⁴C contamination because it involves the introduction of carbon to sample material through the use of laurylamine in the frothing solution. The potential of laurylamine to contaminate in situ ¹⁴C concentrations depends on the carbon source of the compound. With a modern source of carbon, laurylamine has the potential to introduce large quantities, relative to the in situ component, of ¹⁴C to samples. The observed changes in ¹⁴C concentration (seemingly dependent on where quartz was isolated and thus the differing procedures used to isolate quartz) necessitated a systematic study into the potential source and scale of contamination and, if possible, how to efficiently and reliably remove it.

issues and conclude that there must be an unidentified source of ¹⁴C contamination. We suspected that the froth flotation

130 2.0 Systematic Investigation

We isolated quartz from a whole rock sample using five different methods in order to investigate the cause of contamination. The sample selected for this purpose is Caledonian trondhjemite bedrock (Ragnhildstveit et al., 1998) from Utsira, Norway. The sample contains significant feldspar, mica and quartz, making it ideal for use with froth flotation. The ¹⁴C concentration of the sample is irrelevant for the present study; what is important is the ability to observe any potential 135 contamination from the froth flotation procedure. Prior to froth flotation, the sample was crushed, milled, sieved (to isolate the 250 - 500 µm fraction), and magnetically separated. Following magnetic separation, quartz was isolated for aliquot 1 without froth flotation via four days on a shaker table in 5 % HF/HNO3 followed by two days in an ultrasonic bath in 1 % HF/HNO3. Aliquot 1 thus forms a baseline against which the other aliquots are compared. Froth flotation was used with aliquots 2 to 5, which were then etched with different acid mixtures (HF and HF/HNO₃), and varied agitation methods (shaker table and an ultrasonic bath; Table 1). Aliquot 2 spent two days on the shaker table in 5 % HF/HNO3 and two days in the ultrasonic bath in 140 1 % HF/HNO₃, which is the minimum duration of etching that all samples receive in our laboratory. Aliquot 3 also spent two days on the shaker table and two days in the ultrasonic bath but was etched in only HF (5 % on the shaker table and 1 % HF in the ultrasonic bath). Aliquot 3 is essentially our standard procedure but without the inclusion of HNO₃. Aliquots 4 and 5 were not etched on the shaker table and both spent two days in the ultrasonic bath, with the former etched in 1 % HF/HNO₃,

```
145
      and the latter etched in 1 % HF.
```

We extracted carbon from the five quartz aliquots using the Tulane University Carbon Extraction and Graphitization System (TU-CEGS) following the method of Goehring et al. (2019). Quartz is step-heated in the presence of a lithium metaborate (LiBO₂) flux and a high-purity O₂ atmosphere, first at 500 °C for 30 minutes, then at 1100 °C for three hours. The former step is to remove any adsorbed atmospheric CO₂ and combust any carbon derived from handling and dust. Released

150 carbon species from the latter 1100 °C step are oxidised to form CO₂ via secondary hot-quartz-bed oxidation. This is followed by cryogenic collection and purification of the CO₂. Sample yields are measured manometrically (Table 2), and samples are diluted with ¹⁴C-free CO₂. A small aliquot of CO₂ is collected for δ^{13} C analysis, and the remaining CO₂ is graphitised using H₂ reduction over an Fe catalyst (e.g. Southon, 2007). Cathodes containing the graphite were sent to the Woods Hole National Ocean Sciences Accelerator Mass Spectrometry (NOSAMS) to measure ¹⁴C/¹³C isotope ratios (Table 2) relative to NIST





- 155 SRM4990c Oxalic Acid II primary standard. The primary standard was produced in the same graphite reactors used for the unknowns, ensuring full internal normalisation. Stable carbon isotope ratios were measured at the UC-Davis Stable Isotope Facility (Table 2). Repeat measurements of the CRONUS-A interlaboratory comparison standard (Jull et al., 2015) and other samples using the TU-CEGS show that the reproducibility of in situ ¹⁴C measurements is approximately 6 % (Goehring et al., 2019). We therefore present our ¹⁴C concentrations with a conservative 6 % uncertainty as this exceeds the reported analytical
- 160 uncertainty for all of our ¹⁴C measurements. Typical total analytical uncertainties are 1.5 to 2.5 % including the blank correction.

We also measured the carbon isotope ratio of laurylamine to both identify the presence of a modern carbon source for our laurylamine, and to permit a mass balance calculation to quantify the amount of laurylamine left behind after the frothing and etching process. We extracted carbon from 1.9 mg of laurylamine using the TU-CEGS. We used the process regularly used in our laboratory to extract carbon from oxalic acid. This was appropriate given the similarity of the

decomposition temperatures of oxalic acid (~189 °C) and laurylamine (~178 °C). We combusted the sample at 150 °C for ten minutes in ~0.2 MPa (or ~2 atm) of ultra-high purity O₂, after which the temperature was increased to 500 °C to ensure complete combustion. The resulting CO₂ was then cryogenically collected and purified, followed by catalytic reduction via H₂ to graphite. As with the five quartz aliquots, the cathode was sent to NOSAMS to measure the ¹⁴C/¹³C isotope ratio relative to NIST SRM4990c Oxalic Acid II primary standard.

3.0 Results

165

Firstly, the fraction modern (Fm) value of laurylamine is 1.0338 ± 0.0020 , indicative of a modern carbon source.

- Results for the five aliquots are shown in Table 2, with the unit yields and ¹⁴C concentrations also presented in Fig. 2. The total carbon yields for aliquots 1 and 2 are lower than those of aliquots 3 to 5. Aliquots 1 and 2 were isolated without froth flotation
 and with the TUCNL standard procedure (including froth flotation), respectively. Because slightly differing masses of quartz were used for in situ ¹⁴C analysis, a direct comparison can be made using the carbon unit yields (Fig. 2 and Table 2). The unit yield for aliquots 1 and 2 are the same within 1σ error. We observe elevated unit yields for aliquots 3 to 5 relative to those of aliquots 1 and 2.
- As with the unit yields, the ¹⁴C concentration of aliquots 1 and 2 are the same within uncertainties and are distinguishable from the unit yields of aliquots 3 to 5 when using the conservative 6 % uncertainty (Fig. 2). We observe elevated ¹⁴C concentrations for aliquots 3 to 5 relative to those of aliquots 1 and 2, with a particularly high ¹⁴C concentration for aliquot 5 (Fig. 2B). Figure 2 shows that the higher unit yields correspond with higher measured ¹⁴C concentrations. With aliquot 5, a small increase in unit yield results in a disproportionately high ¹⁴C concentration that dwarfs those of aliquots 1 to 4.





185 4.0 Discussion

The modern carbon source, identified with the measured Fm of our laurylamine, shows that laurylamine is not ¹⁴C dead and thus has the potential to contaminate samples with respect to ¹⁴C. The measured ¹⁴C/¹²C ratio for laurylamine is 1.19 x 10⁻¹². This means that, for example, 20 μg contains ~9.3 x 10⁵ atoms of ¹⁴C. The elevated unit yields and ¹⁴C concentrations of aliquots 3 to 5 relative to those of aliquots 1 and 2 indicate that the former are contaminated with total carbon and, of particular importance, ¹⁴C. The unit yields and ¹⁴C concentrations of aliquots 1 and 2 are indistinguishable from one another, which indicates that our standard procedure for quartz isolation (Aliquot 2) removes carbon introduced by laurylamine. Differing quartz isolation procedures used at other laboratories may therefore explain why quartz isolated from the same samples at Tulane and elsewhere produced vastly different ¹⁴C concentrations and unit yields (Sect. 1.2).

We use the excess measured ¹⁴C atoms in aliquots 3 to 5 (the total ¹⁴C atoms for each aliquot in excess of the average of those of aliquots 1 and 2) with the measured ¹⁴C/¹²C ratio for laurylamine to calculate the corresponding mass of residual carbon and laurylamine, per gram of quartz, that was collected with the in situ ¹⁴C component. We assume aliquots 1 and 2 were not contaminated with ¹⁴C and assume that the excess ¹⁴C is sourced solely from laurylamine. To calculate the mass of contaminant (*M_{contam}*) carbon or laurylamine we follow

200
$$M_{contam} = {}^{14}C_{excess} \left(\frac{{}^{14}C}{{}^{12}C}\right)_{LA} \frac{M}{A}$$

where ${}^{14}C_{excess}$ is the measured number of excess ${}^{14}C$ atoms, $({}^{14}C/{}^{12}C)_{LA}$ is the measured ratio for laurylamine (1.19 x 10⁻¹²), M is the molar mass of carbon or molecular mass of laurylamine, and A Avogadro's Number. This calculation is an estimate as it does not take into account the ca. 1.1 % ${}^{13}C$ in laurylamine. The excess ${}^{14}C$ accounts for an estimated 0.06, 0.03 and 0.74 µg carbon g^{-1} quartz, and 0.08, 0.04, 0.95 µg g^{-1} of laurylamine for aliquots 3, 4 and 5, respectively (Table 2). We can use the same method for the samples that produced the initial anomalous measurements shown in Fig. 1. To do so we assume that the final measurement made for each sample is free from laurylamine contamination. For the samples presented in Fig. 1, the excess ${}^{14}C$ concentrations range from 1.38 x 10⁵ to 3.23 x 10⁵ at g^{-1} . The associated residual carbon ranges from 2.32 to 5.42 µg g^{-1} , and the residual laurylamine ranges from 2.98 to 6.96 µg g^{-1} , both per gram of quartz. We speculate that the latter

210 residual carbon and laurylamine estimates, an order of magnitude greater than those presented in this study, may be an artefact of the differing froth flotation and etching procedures used at the laboratories at which the quartz was isolated. Contributing factors could include, but are not limited to, a greater amount of laurylamine used in the quartz separation process, the concentration at which the laurylamine comes into contact with sample material (dilute or undilute), the acids used in the etching procedure, and the duration of acid etching.

The elevated carbon yields and ¹⁴C concentrations of aliquots 3 to 5 relative to those of aliquots 1 and 2 suggest two things. Firstly, it is apparent that HNO₃ is needed to remove laurylamine-derived carbon, both total carbon and ¹⁴C,





contamination from quartz. The importance of HNO3 is demonstrated by the higher unit yields and ¹⁴C concentrations of aliquots that were etched with only HF compared to aliquots that had the same quartz isolation method and duration of etching but were etched with a combination of HF and HNO3 (aliquot 3 vs 2 and aliquot 5 vs 4). We hypothesise that, as an oxidiser, 220 HNO₃ is key in the decomposition of laurylamine. Before carbon is extracted from aliquots, quartz is etched for 30 minutes in 50 % HNO₃ (Lifton et al., 2001; Goehring et al., 2019). This is important to note because it is apparent that this additional etch with strong HNO₃ is not sufficient alone to remove contaminant ¹⁴C and highlights the importance of HF as well as HNO₃ in the etching procedure and their role in the removal of contamination. We speculate that dissolution of quartz using HF helps to release laurylamine stored within microfractures of quartz grains (elaborated further below). Secondly, two days in an 225 ultrasonic bath, regardless of whether HF or HF/HNO3 is used for etching, appears to be insufficient to remove contaminant laurylamine. Aliquots 4 and 5, which were not etched on the shaker table and spent a total of two days etching in an ultrasonic bath, both appeared visually pure and thus looked ready for in situ ¹⁴C analysis without the context of potential laurylamine contamination. A standard procedure to etch samples for subsequent in situ¹⁴C analysis until they are visually pure is evidently not sufficient when froth flotation has been used.

230

The fact that the ¹⁴C concentration increase from laurylamine is of the same order of magnitude as that of typical in situ ¹⁴C measurements is of great concern and highlights the need for a standardised procedure to eliminate contamination from the quartz isolation process. Carbon introduced by laurylamine is evidently persisting through a 500 °C bake, the first stage of extracting ¹⁴C from quartz with the TU-CEGS and in other in situ ¹⁴C laboratories (e.g., Hippe et al., 2013; Lifton et al., 2015; Goehring et al., 2019; Lamp et al., 2019). The 500 °C bake was previously shown to remove contaminant ¹⁴C (Lifton

- 235 et al, 2001), though this was from sample handling and the atmosphere and pre-dates the implementation of froth floatation for quartz separation. We suspect that the observed contamination is sourced from laurylamine residing within microfractures of quartz grains, which may explain why the contamination is able to persist through the 500 °C bake and possibly accounts for differences in the degrees of contamination between previously analysed samples and those as part of this study due to differences in the quartz grain characteristics. Further etching, both in our initial measurements (Sect. 1.2) and when comparing
- aliquots 2 and 3 with aliquots 4 and 5, lowers unit yields and ¹⁴C concentrations. The longer duration in acid may indicate that 240 the HF is opening up microfractures and freeing contamination to be removed, highlighting the importance of HF in the removal of contamination, though this would be difficult to test.

5.0 Conclusion

We found that laurylamine is manufactured with a modern carbon source and thus introduces modern ¹⁴C to sample 245 material during froth flotation. We have shown through a systematic study that contaminant ¹⁴C from laurylamine persists through sample etching and is collected with in situ ¹⁴C if etching is not rigorous enough. Nitric acid, combined with hydrofluoric acid, is required to effectively remove contaminant ¹⁴C, which is shown by the elevated ¹⁴C concentrations of guartz separates isolated without nitric acid relative to those extracted with nitric acid. We have outlined a reliable method for





ensuring no contaminant ¹⁴C from laurylamine remains with quartz following froth flotation and etching. In short, two 24-hour
 periods on a shaker table with 5 % HF/HNO₃, followed by two 24-hour periods in an ultrasonic bath in 1 % HN/HNO₃, is sufficient to produce in situ ¹⁴C concentrations indistinguishable from a sample for which quartz was isolated without froth flotation. Ultimately, froth flotation should be used with caution and the sample etching procedure outlined above should be used at an absolute minimum.

Data availability

255 All data are included with this manuscript.

Author contributions

The study was conceived by BG and KN. Sample collection, preparation and analysis undertaken by KN and BG. Manuscript was written by KN and BG.

Competing interests

260 The authors declare that they have no conflict of interest.

Acknowledgements

KN and BG would like to acknowledge NOSAMS for their outstanding measurements. We would also like to acknowledge the fruitful discussions with Trevor Hillenbrand, John Stone, and Greg Balco that greatly helped to develop the study. Sample material was collected with the help of Jan Mangerud and John-Inge Svendsen.

265 **References**

Balco, G., Todd, C., Goehring, B.M., Moening-Swanson, I., Nichols, K.: Glacial geology and cosmogenic-nuclide exposure ages from the Tucker Glacier - Whitehall Glacier confluence, northern Victoria Land, Antarctica, Am. J. Sci., 319(April), 255–286, doi:10.2475/04.2019.01, 2019.

Balco, G., Todd, C., Huybers, K., Campbell, S., Vermeulen, M., Hegland, M., Goehring, B. M. and Hillebrand, T. R.:

270 Cosmogenic-nuclide exposure ages from the Pensacola Mountains adjacent to the foundation ice stream, Antarctica, Am. J. Sci., 316(6), 542–577, doi:10.2475/06.2016.02, 2016.

Goehring, B. M., Wilson, J. and Nichols, K.: A fully automated system for the extraction of in situ cosmogenic carbon-14 in the Tulane University cosmogenic nuclide laboratory, Nucl. Instruments Methods Phys. Res. Sect. B Beam Interact. with Mater. Atoms, doi:10.1016/j.nimb.2019.02.006, 2019.

Herber, L. J.: Separation of feldspar from quartz by flotation, Am. Mineral., 54, 1212–1215, doi:10.4144/rpsj1954.25.192, 1969.

Hippe, K., Kober, F., Wacker, L., Fahrni, S. M., Ivy-Ochs, S., Akçar, N., Schlüchter, C. and Wieler, R.: An update on in situ





cosmogenic 14C analysis at ETH Zürich, Nucl. Instruments Methods Phys. Res. Sect. B Beam Interact. with Mater. Atoms, 294, 81–86, doi:10.1016/j.nimb.2012.06.020, 2013.

- Jull, A. J. T., Scott, E. M. and Bierman, P.: The CRONUS-Earth inter-comparison for cosmogenic isotope analysis, Quat. Geochronol., 26(1), 3–10, doi:10.1016/j.quageo.2013.09.003, 2015.
 Kohl, C. P. and Nishiizumi, K.: Chemical isolation of quartz for measurement of in-situ -produced cosmogenic nuclides, Geochim. Cosmochim. Acta, 56(9), 3583–3587, doi:10.1016/0016-7037(92)90401-4, 1992.
 Lamp, J. L., Young, N. E., Koffman, T., Schimmelpfennig, I., Tuna, T., Bard, E. and Schaefer, J. M.: Update on the
- cosmogenic in situ 14C laboratory at the Lamont-Doherty Earth Observatory, Nucl. Instruments Methods Phys. Res. Sect. B Beam Interact. with Mater. Atoms, (April), 1–6, doi:10.1016/j.nimb.2019.05.064, 2019.
 Lifton, N., Goehring, B., Wilson, J., Kubley, T. and Caffee, M.: Progress in automated extraction and purification of in situ 14C from quartz: Results from the Purdue in situ 14C laboratory, Nucl. Instruments Methods Phys. Res. Sect. B Beam Interact. with Mater. Atoms, 361, 381–386, doi:10.1016/j.nimb.2015.03.028, 2015.
- 290 Lifton, N. A., Jull, A. J. T. and Quade, J.: A new extraction technique and production rate estimate for in situ cosmognic 14C in quartz, Geochim. Cosmochim. Acta, 65(12), 1953–1969, doi:10.1016/S0016-7037(01)00566-X, 2001. Ragnhildstveit, E., Naterstad, J., Jorde, K., Egeland, B.: Beskrivelse Til Geologisk Kart Over Norge - 1: 250.000 Haugesund, Norges Geol. Undersøkelse., 1998.
- Southon, J.: Graphite reactor memory Where is it from and how to minimize it?, Nucl. Instruments Methods Phys. Res.
- 295 Sect. B Beam Interact. with Mater. Atoms, 259(1), 288–292, doi:10.1016/j.nimb.2007.01.251, 2007.

Figures

Figure 1: Elevation versus in situ ¹⁴C concentration of samples that initially yielded anomalously high ¹⁴C concentrations. Measurements at the same elevation are from the same sample. For each sample, the highest in situ ¹⁴C concentration is sourced

- 300 from the first measurement. For the two samples measured three times at 510 and 875 m asl, the second measurement was made following additional etching and yields the intermediate in situ ¹⁴C concentration. For every sample, the final lowest concentration measurement was made from quartz isolated from whole rock at Tulane. Error bars reflect a long-term 6 % uncertainty. Some error bars are smaller than their respective data points. Thick grey line and grey shading are the saturation concentration and associated error envelope.
- 305 Figure 2: A: Unit yield for the five quartz separates. Aliquot numbers refer to those in Table 1. B and C: C-14 concentration for the five quartz aliquots; error bars reflect a long-term 6 % uncertainty. Both B. and C. present the same ¹⁴C concentrations except that C. has a smaller axis range to highlight differences between the first four aliquots. See Table 1 for the different quartz isolation procedures used. For reference, all aliquots other than Aliquot 1 were subject to froth flotation. Aliquot 2 was processed using the TUCNL standard procedure.

310 Tables

Table 1: Aliquot information and quartz isolation procedures. Whilst on the shaker table, samples were etched in 5 % HF/HNO₃ or 5 % HF. Whilst in the ultrasonic bath, samples were etched in 1 % HF/HNO₃ or 1 % HF.

Table 2: In situ ¹⁴C analytical data. Effective blank is representative of the blank during the running of the samples presented. SeeSect. 2 for rationale behind the use of the 6 % error for the ¹⁴C concentrations. The mass of residual carbon and laurylamine are315calculated using the measured ¹⁴C/¹²C ratio of laurylamine (see Sect. 4).







Fig. 1









325





| Sample ID | Aliquot Number | Days on shaker table | Days in ultrasonic bath | Total | Notes |
|--------------------------|-------------------|----------------------------|-------------------------------|-------|--|
| 16-UT-004-QUA-NOFROTH | 1 | 4 | 2 | 6 | No froth flotation, etched with HF/HNO ₃ |
| 16-UT-004-QUA-NORM | 2 | 2 | 2 | 4 | Froth flotation, etched with HF/HNO ₃ in shaker table and ultrasonic bath |
| 16-UT-004-QUA-HFONLY | 3 | 2 | 2 | 4 | Froth flotation, etched with HF in shaker table and ultrasonic bath |
| 16-UT-004-QUA-NOSTABLE1 | 4 | 0 | 2 | 2 | Froth flotation, etched with HF/HNO3 in ultrasonic bath |
| 16-UT-004-QUA-NOSTABLE12 | 5 | 0 | 2 | 2 | Froth flotation, etched with HF in ultrasonic bath |

Table 1.

335

340





| Aliquot Number | TUCNL | AMS Lab | AMS ID | Quartz weight | C yield | ± | Unit Yield | Diluted Gas Mass | ±355 |
|-------------------|-------|---------|---------------|------------------|------------|------|--------------------|------------------------|------------|
| | | | | (g) | (µg) | (µg) | $(\mu g \ g^{-1})$ | (µg) | (µg) |
| 1 | 309 | NOSAMS | OS- 141782 | 5.196 | 6.9 | 0.1 | 1.3 | 85.5 | 1.1 |
| 2 | 307 | NOSAMS | OS- 141784 | 5.122 | 6.5 | 0.1 | 1.3 | 100.3 | 1.3 |
| 3 | 308 | NOSAMS | OS- 141785 | 5.104 | 8.3 | 0.1 | 1.6 | 88.0 | 360 1.1 |
| 4 | 310 | NOSAMS | OS- 141786 | 5.134 | 8.4 | 0.1 | 1.6 | 83.7 | 1.1 |
| 5 | 311 | NOSAMS | OS- 141788 | 5.080 | 9.1 | 0.1 | 1.8 | 76.7 | 1.0 |

365

| ¹⁴ C/ ¹³ C corrected | 1σ error | $\delta^{13}C$ | Cor. \pm | ¹⁴ C/C total | 1σ error | Total ¹⁴ C atoms | 1σ error | ¹⁴ C conc. |
|---|-----------------|----------------|------------|----------------------------|-----------------|--------------------------------|-----------------|-----------------------|
| | | (‰) | | | | blank corrected | (at) | (at.g ⁻¹) |
| 8.47E-12 | 8.56E-14 | -4.98 | 0.5 | 9.29E-14 | 9.39E-16 | 3.34E+05 | 6.85E+03 | 6.42E+04 |
| 6.88E-12 | 8.23E-14 | -4.54 | 0.5 | 7.55E-14 | 9.04E-16 | 3.15E+05 | 6.85E+03 | 6.15E+04 |
| 9.88E-12 | 9.58E-14 | -5.20 | 0.5 | 1.08E-13 | 1.05E-15 | 4.13E+05 | 6.85E+03 | 8.10E+04 |
| 9.57E-12 | 9.27E-14 | -5.13 | 0.5 | 1.05E-13 | 1.02E-15 | 3.76E+05 | 6.85E+03 | 7.32E+04 |
| 3.64E-11 | 1.90E-13 | -5.78 | 0.5 | 3.99E-13 | 2.08E-15 | 1.47E+06 | 6.85E+03 | 2.89E+05 |

| 1σ error | 6% Error | Effective blank | lσ error | Effective Blank as % | Residual C | Residual C12H27N |
|-----------------------|-----------------------|-----------------|----------|---------------------------------------|--------------------|---------------------|
| (at.g ⁻¹) | (at.g ⁻¹) | (at) | (at) | of total ¹⁴ C At Sample | $(\mu g \ g^{-1})$ | $(\mu g \ g^{-1})$ |
| 1.64E+03 | 3.85E+03 | 6.47E+04 | 6.85E+03 | 16.25 | - | - |
| 1.64E+03 | 3.69E+03 | 6.47E+04 | 6.85E+03 | 17.04 | - | - |
| 1.80E+03 | 4.86E+03 | 6.47E+04 | 6.85E+03 | 13.53 | 0.06 | 0.08 |
| 1.73E+03 | 4.39E+03 | 6.47E+04 | 6.85E+03 | 16.25 | 0.03 | 0.04 |
| 4.10E+03 | 1.74E+04 | 6.47E+04 | 6.85E+03 | 14.69 | 0.74 | 0.95 |

370

Table 2. (split into three sections)